



eed a lightweight folding wheelchair with simple adjustability and reliable performance at an economical price? The Invacare® Patriot™ is the model of choice. The classic manual rehab wheelchair is a clean design offering valuable features including: a lightweight, aluminum frame; adjustable seat to floor heights; durable center mount front riggings; and adjustable angle head tubes all in a compact folding design.



Simple, triangular axle bracket offers 8 different axle positions – 4 axle positions in font of back canes and 4 axle positions behind back canes.

Highlights

- Lightweight, aluminum frame
- Adjustable seat-to-floor height and rear axle positions
- · Adjustable angle head tubes
- Durable center mount front riggings
- Extra smooth precision head tube bearings
- Compact for transport (11.75" wide folded with wheels)
- Adjustable-height back available in I" increments, 15"-19"
- Multiple arm, footplate and tire options
- 14 color options
- Quick-release axle option
- · Seat pouch standard



Precision head tube bearings allow casters to turn smoothly.

Specifications

| Overall | | | | |
|---------------------|-------------------------|--------------|------------|--|
| Width open | 7.50" plus seat width* | | | |
| Width closed | I I.75"* | | | |
| Depth (w/riggings) | 43" or 45" | | | |
| Seat widths | 14", 16", 17", 18", 20" | | | |
| Seat depths | 16", 18" | | | |
| Front Seat-to-Floor | | Front Caster | Rear Wheel | |

| it seat-to-rioor | | | Front Caster | Real vvilleel | |
|-------------------|---------------|---------------------------------------|--------------|-------------------------|-------------|
| D \A# 16: | C | • | Mounting | Mounting Po | |
| Rear Wheel Size | Seat-to-Floor | Caster | Position | Axle Mounti | ng Bracket |
| 20" | 15.50" | 6" | top | middle | upper |
| | 16.50" | 6" | † | refer to service manual | |
| | 17.50" | 6" | bottom | bottom | lower |
| | 17.50" | 8" | top | bottom | lower |
| | 18.50" | 8" | t | refer to serv | vice manual |
| 22" | 15.50" | 6" | † | refer to serv | vice manual |
| | 16.50" | 6" | † | refer to service manual | |
| | 17.50" | 6" | bottom | bottom | lower |
| | 17.50" | 8" | top | bottom | lower |
| | 18.50" | 8" | t | refer to serv | vice manual |
| | 19.50" | 8" | t | refer to service manual | |
| 24" | 16.50" | 6" | † | refer to service manual | |
| | 17.50" | 6" | bottom | middle | upper |
| | 17.50" | 8" | top | middle | upper |
| | 18.50" | 8" | t | refer to serv | vice manual |
| | 19.50" | 8" | bottom | bottom | lower |
| nee-to-heel** | Min. | Max. (range adjusts in 1" increments) | | | |
| 70 degree tapered | 14.75" | 18.75" | | | |
| 70 degree | 12.50" | 17.50" | | | |
| 60 degree | 13.75" | 17.50" | | | |
| Elevating | 14.50" | 20.25" | | | |
| | | | | | |

Arm Types

Conventional dual-point

Desk or Full Length

Fixed or Adjustable Height

Space-Saver dual-point

Desk or Full Length

| Fixed or Adjustable Height | | | |
|----------------------------|-----------------------------------|--|--|
| Back heights | 15"-19" | | |
| Weight capacity | 250 lb. | | |
| Weight | 29 lb. | | |
| Limited warranty | Lifetime on frame and crossbraces | | |

 $^{{}^*\!}M\!$ easurements taken from outside of handrim to outside of handrim with 0 degrees camber

 $\dagger Not$ available as a factory setting



Invacare Corporation www.invacare.com

USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

Canada

570 Matheson Blvd. E. Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

Invacare is a registered trademark and Yes, you can., the blue medallion design and Patriot are trademarks of Invacare Corporation.

Specifications are subject to change without notification.

© 2003 Invacare Corporation

Form No. 03-051



"Invacare says 'Yes, you can",' and I couldn't have said it betten myself."

Servold talmer

^{**}Measurement taken from top of seat rail to back of footrest